



piccolo^{xpress} chemistry analyzer

FAST, EASY, ACCURATE ON-SITE CHEMISTRY TESTING

URGENT CARE MEANS BEING PREPARED TO DIAGNOSE AND TREAT A WIDE VARIETY OF CONDITIONS SWIFTLY, ACCURATELY AND RELIABLY.



EXCEPTIONALLY VERSATILE

Wide range of **CLIA waived** panels to increase your on-site diagnostic capabilities

- **Comprehensive Metabolic Panel**
ALB, ALP, ALT, AST, BUN, Ca, Cl⁻, CRE, GLU, K⁺, Na⁺, TBIL, tCO₂, TP
- **Basic Metabolic Panel**
BUN, Ca, Cl⁻, CRE, GLU, K⁺, Na⁺, tCO₂
- **Lipid Panel**
CHOL, CHOL/HDL*, HDL, LDL*, TRIG, VLDL*
- **Lipid Panel Plus**
ALT, AST, CHOL, CHOL/HDL*, GLU, HDL, LDL*, TRIG, VLDL*
- **Liver Panel Plus**
ALB, ALP, ALT, AMY, AST, GGT, TBIL, TP
- **General Chemistry 6**
ALT, AST, BUN, CRE, GGT, GLU
- **General Chemistry 13**
ALB, ALP, ALT, AMY, AST, BUN, Ca, CRE, GGT, GLU, TBIL, TP, UA
- **Electrolyte Panel**
Cl⁻, K⁺, Na⁺, tCO₂
- **Kidney Check**
BUN, CRE
- **Renal Function Panel**
ALB, BUN, Ca, Cl⁻, CRE, GLU, K⁺, Na⁺, PHOS, tCO₂
- **MetLyte 8 Panel**
BUN, CK, Cl⁻, CRE, GLU, K⁺, Na⁺, tCO₂

- **Hepatic Function Panel**
ALB, ALP, ALT, AST, DBIL, TBIL, TP
- **Basic Metabolic Panel Plus**
BUN, Ca, Cl⁻, CRE, GLU, K⁺, Na⁺, tCO₂, Mg, Lactate Dehydrogenase
- **MetLyte Plus CRP**
BUN, CK, Cl⁻, CRE, GLU, K⁺, Na⁺, tCO₂, CRP

○ MAXIMIZE PATIENT THROUGHPUT

Fast, lab-accurate and CLIA waived general chemistry results

- Comprehensive health screens within the span of each patient visit
- 14 complete panels, 11 of which are CLIA waived
- Diagnose nearly any patient condition and either treat or triage
- Confident clinical diagnostic information in just 12 minutes
- Maximize the value of every patient visit

○ PROVIDE COMPREHENSIVE CARE TO THE PATIENTS YOU SERVE

Be prepared to diagnose nearly any patient condition that is presented

- Quickly screen in or rule out clinical indications on-site in minutes
- Diagnose and triage complex metabolic, kidney, liver and cardiovascular disorders
- Serve a greater audience of afflictions, including hypertension and diabetes
- Treat abdominal pain with accurate diagnoses from accurate chemistry results
- Monitor the effects of prescribed medications – everything from statins to diabetes therapy
- Weakness, fatigue and dizziness are a sign of electrolyte deficiencies
- Expand your offering with up to 29 different analytical chemistry results

CLIA WAIVED

*Calculated



3240 Whipple Road, Union City, CA 94587 | 1-800-822-2947 | www.abaxis.com

©2010 Abaxis, Inc., Piccolo Xpress™ is a trademark of Abaxis, Inc. 888-3223 Rev.A

PICCOLO XPRESS PANEL SPECIFICATIONS

| | <ul style="list-style-type: none"> • Comprehensive Metabolic Panel • Basic Metabolic Panel • Lipid Panel • Lipid Panel Plus • Liver Panel Plus • General Chemistry 6 • General Chemistry 13 • Electrolyte Panel • Kidney Check • Renal Function Panel • MetLyte 8 Panel • Hepatic Function Panel • Basic Metabolic Panel Plus • MetLyte Plus CRP | | | | | | | | | | | | | |
|------------------|--|---|---|---|---|---|---|---|---|---|--------------------|---|---|---|
| | CLIA WAIVED | | | | | | | | | | MODERATELY COMPLEX | | | |
| ALB | ● | | | | ● | | ● | | | ● | | ● | | |
| ALP | ● | | | | ● | | ● | | | | | ● | | |
| ALT | ● | | | ● | ● | ● | ● | | | | | ● | | |
| AMY | | | | ● | | ● | | | | | | | | |
| AST | ● | | | ● | ● | ● | | | | | | ● | | |
| DBIL | | | | | | | | | | | | ● | | |
| BUN | ● | ● | | | | ● | ● | | ● | ● | ● | | ● | ● |
| Ca | ● | ● | | | | | ● | | ● | | | | ● | |
| CK | | | | | | | | | | | ● | | | ● |
| Cl ⁻ | ● | ● | | | | | | ● | ● | ● | | | ● | ● |
| CRE | ● | ● | | | | ● | ● | | ● | ● | | | ● | ● |
| CRP | | | | | | | | | | | | | | ● |
| GGT | | | | ● | ● | ● | | | | | | | | |
| GLU | ● | ● | | ● | ● | ● | | | ● | ● | | | ● | ● |
| K ⁺ | ● | ● | | | | | | ● | ● | ● | | | ● | ● |
| LDH | | | | | | | | | | | | | ● | |
| Mg | | | | | | | | | | | | | ● | |
| Na ⁺ | ● | ● | | | | | | ● | ● | ● | | | ● | ● |
| PHOS | | | | | | | | | ● | | | | | |
| TBIL | ● | | | ● | | ● | | | | | | ● | | |
| tCO ₂ | ● | ● | | | | | | ● | ● | ● | | | ● | ● |
| TP | ● | | | ● | | ● | | | | | | ● | | |
| UA | | | | | | ● | | | | | | | | |
| CHOL | | | ● | ● | | | | | | | | | | |
| CHOL/HDL* | | | ● | ● | | | | | | | | | | |
| HDL | | | ● | ● | | | | | | | | | | |
| LDL* | | | ● | ● | | | | | | | | | | |
| TRIG | | | ● | ● | | | | | | | | | | |
| VLDL* | | | ● | ● | | | | | | | | | | |

*Calculated value